


# Product information

## KYPHOKIT 11G

Fast access (no k-wires needed)



## General information

---

- **Cost efficient:** individual set configuration possible
- Save one step as opposed to conventional kyphoplasty systems („**fast access**“): balloon can be inserted directly through the needle (no k-wire)

## Features/Design

---

- Access: 11G (or 8G -> see also KYPHOKIT 8G)
- Familiar handling: KYPHOKIT 11G conforms with other market based systems for kyphoplasty

# Order information

## KYPHOKIT 11G



### KYPHOKIT 11G Sets

Product name	Ref-No.	Description
KYPHOKIT 11G 15	KK11-101V-15J	KYPHOKIT 11G with 15 mm balloon catheter (400PSI/27ATM) also included: 1x access needle, 1x syringe pump, 1x drill, 3x bone filler, 4x cement syringe
KYPHOKIT 11G 20	KK11-101V-15J	KYPHOKIT 11G with 20 mm balloon catheter (400PSI/27ATM) also included: 1x access needle, 1x syringe pump, 1x drill, 3x bone filler, 4x cement syringe

### KYPHOKIT 11G Individual Components

Product name	Ref-No.	Description
balloon catheter 15mm	KBAL15J	balloon catheter 15mm (400PSI/27ATM)
balloon catheter 20mm	KBAL20J	balloon catheter 20mm (400PSI/27ATM)
syringe pump	K-INFA	analogue syringe pump
drill	KK11-DRI	drill 11G
bone filler	KK11-FIL	bone filler 11G, charge 0.7 ml, Luer Lock
access needle beveled tip	KK11-WB	fast access needle with beveled tip 11G, Luer Lock
access needle trocar tip	KK11-WTR	fast access needle with trocar tip 11G working length 120mm + 3-4cm handle, Luer Lock